

## Tapones de corcho



- Para varios usos
- Diferentes diámetros disponibles

Ref.	Ø inferior (mm)	Ø superior (mm)	Altura (mm)	Unid./caja	€/caja
764015	3	6	16	100	NC -
764016	5	8	16	100	NC -
764017B	6	9	17	100	NC -
764019	10	13	22	100	NC -
764020	11	14	22	100	NC -
764021	13	16	22	100	NC -
764022	15	18	22	100	NC -
764023	17	20	26	50	NC -
764024	18	22	26	50	NC -
764025	26	30	27	50	NC -
764026	32	36	27	50	NC -
764027	36	40	27	20	NC -

## Gradillas de acero inoxidable autoclavables con revestimiento RILSAN

Ref.	L x An x Al (mm)	Para tubos de ø (mm)	Número de tubos	€
805108*	207 x 103 x 25	11	5 x 10	NC -
844038	72 x 55 x 55	13	4 x 3	NC -
844068	108 x 72 x 55	13	4 x 6	NC -
844128	206 x 72 x 55	13	4 x 12	NC -
840108	172 x 172 x 55	13	10 x 10	NC -
854038	100 x 78 x 65	16 - 18 - 20	4 x 3	NC -
864038	100 x 78 x 100	16 - 18 - 20	4 x 3	NC -
864068	150 x 100 x 100	16 - 18 - 20	25	NC -
864128	296 x 100 x 100	16 - 18 - 20	4 x 12	NC -
874128	348 x 122 x 100	22 - 25	4 x 12	NC -
843288	280 x 110 x 90	30	3 x 8	NC -

\* Gradilla con asas laterales.

## Gradillas de acero inoxidable



Ref.	L x An x Al (mm)	Para tubos de ø (mm)	Número de tubos	€
844061	108 x 72 x 55	13	4 x 6	NC -
844121	206 x 72 x 55	13	4 x 12	NC -
854061	150 x 100 x 65	16 - 18 - 20	4 x 6	NC -
864061	150 x 100 x 100	16 - 18 - 20	4 x 6	NC -
864121	296 x 100 x 100	16 - 18 - 20	4 x 12	NC -
874061	176 x 122 x 100	22 - 25	4 x 6	NC -
874121	348 x 122 x 100	22 - 25	4 x 12	NC -
841121	205 x 20 x 50	tubo de hemólisis	1 x 12	NC -
843201	291 x 54 x 50	45	1 x 6	NC -

## Media gradilla de acero y epoxi



Gradilla pequeña de acero y epoxi ideal para los espacios pequeños: bajo una campana, al baño maría, etc.  
Disponible en 3 colores: verde, azul o naranja.

Para tubos Ø (mm)	Capacidad tubos	L x An x Al (mm)	Verde	Azul	Naranja	€
10 a 13	36 (6 x 6)	110 x 130 x 62	-	086609	-	NC -
15 a 16	24 (6 x 4)	121 x 89 x 62	086607	086610	086613	NC -
18 a 20	20 (5 x 4)	125 x 108 x 62	086608	086611	086614	NC -